

Rad17 (phospho Ser645) antibody [EP1519Y]

Cat. No. GTX61991

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at a dilution of 1:2500 - 5000

Optimal working dilution for a specific application should be determined by the investigator.

Calculated MW	77 kDa. (Note)
Product Note	The antibody only detects Rad17 phosphorylated on Serine 645.


Properties

Form	Liquid
Buffer	50mM Tris-Glycine, 150mM NaCl, 0.05% BSA, 40% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A phospho-specific peptide corresponding to residues surrounding Serine 645 of human Rad17 was used as immunogen.
Purification	Tissue culture supernatant
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. RabMAb® technology is covered by the following U.S. Patents, No. 5,675,063 and/or 7,429,488.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX61991 WB Image

A. Western blot analysis on HeLa cell lysates using anti-Rad17 phospho (pS645) RabMAb (cat. #GTX61991) 1:5000 dilution. Cells were either (A) untreated (B) treated with Lambda Phosphatase.



For full product information, images and publications, please visit our [website](https://www.genetex.com).