

RIAM antibody [EPR2806]

Cat. No. GTX62255

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at a dilution of 1:2000 - 10,000

For IHC: Use at a dilution of 1:100 - 250

For ICC: Use at a dilution of 1:100 - 250

For IP: Use at a dilution of 1:50

Optimal working dilution for a specific application should be determined by the investigator.

Calculated MW 73 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	50mM Tris-Glycine, 150mM NaCl, 0.05% BSA, 40% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A synthetic peptide corresponding to residues in human RIAM was used as an immunogen.
Purification	Tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx62255 WB Image

Western blot analysis on (A) Human Lymph node (B) Fetal brain (C) human spleen (D) Jurkat and (E) HL-60 + retinoic acid cell lysates using anti-RIAM RabMAb (cat. #GTx62255).



GTx62255 IHC-P Image

Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-RIAM RabMAb (cat. #GTx62255-1).



For full product information, images and publications, please visit our [website](https://www.genetex.com).