

ING4 antibody [EP3804]

Cat. No. GTX62655

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human, Mouse, Rat

References (2)

Package

100 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at a dilution of 1:1,000 - 10,000

For IP: Use at a dilution of 1:10 - 100

Optimal working dilution for a specific application should be determined by the investigator.


Calculated MW 29 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	50mM Tris-Glycine, 150mM NaCl, 0.05% BSA, 40% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A synthetic peptide corresponding to residues in human ING4 was used as an immunogen.
Purification	Tissue culture supernatant
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. RabMAb® technology is covered by the following U.S. Patents, No. 5,675,063 and/or 7,429,488.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx62655 WB Image

Western blot analysis on (A) 293T (B) T47-D and (C) Jurkat cell lysates using anti-ING4 RabMAb (cat. #GTx62655).


GTx62655 WB Image

Western blot analysis on (A) C6 (B) Raw264.7 and (C) PC-12 cell lysates using anti-ING4 RabMAb (cat. #GTx62655).



For full product information, images and publications, please visit our [website](https://www.genetex.com).