

CD80 antibody [EPR1157(2)]

Cat. No. GTX62858

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µl

Applications

Application Note

Is unsuitable for ICC/IF or IHC-P. WB: 1/1000 - 1/10000. Detects a band of approximately 60 kDa.

Calculated MW 33 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	Buffer, 0.05% BSA, 50% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to a region within human CD80.
Purification	Protein A purified
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. RabMAb® technology is covered by the following U.S. Patents, No. 5,675,063 and/or 7,429,488.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx62858 WB Image

WB analysis of Raji lysate (20 µg) using CD80 antibody [EPR1157(2)] at a dilution of 1:2500.



For full product information, images and publications, please visit our [website](#).