

HSV tag antibody [GT212]

Cat. No. GTX628790

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF
Reactivity	Species independent

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen used to generate this antibody corresponds to Hsv Tag
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX628790 WB Image


Hsv Tag antibody [GT212] detects Hsv protein by Western blot analysis.

A. 30 µg 293T whole cell lysate/extract

B. 30 µg whole cell lysate/extract of HSV-tagged protein expressing 293T cells

12 % SDS-PAGE

Hsv Tag antibody [GT212] (GTX628790) dilution: 1:5000



GTX628790 ICC/IF Image

HSV-tag antibody detects HSV-recombinant protein by immunofluorescent analysis. Sample: EGFP/HSV-tag transfected HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Red: HSV-recombinant protein stained by HSV-tag antibody (GTX628790) diluted at 1:2000. Green: GFP is expressed in the transfected cell. Blue: Hoechst 33342 staining.



For full product information, images and publications, please visit our [website](https://www.genetex.com).