

ZNF703 antibody [GT264]

Cat. No. GTX628896

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 58 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.21 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human ZNF703. The exact sequence is proprietary.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX628896 ICC/IF Image

ZNF703 antibody [GT264] detects ZNF703 protein at nucleus by immunofluorescent analysis.
 Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.
 Green: ZNF703 stained by ZNF703 antibody [GT264] (GTX628896) diluted at 1:2000.
 Blue: Hoechst 33342 staining.
 Scale bar = 10 μm.



GTX628896 WB Image

ZNF703 antibody detects ZNF703 protein by western blot analysis. Non-transfected (-) and ZNF703-transfected (+,) 293T whole cell extracts (30 μg) were separated by 10% SDS-PAGE, and the membrane was blotted with ZNF703 antibody (GTX628896) at a dilution of 1:5000.



For full product information, images and publications, please visit our [website](https://www.genetex.com).