

Factor X antibody [GT268]

Cat. No. GTX629480

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:5000-1:20000
ICC/IF	1:100-1:1000

Not tested in other applications.

Calculated MW 55 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human Factor X. The exact sequence is proprietary.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx629480 ICC/IF Image

Factor X antibody [GT268] detects Factor X protein by immunofluorescent analysis.
 Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.
 Green: Factor X stained by Factor X antibody [GT268] (GTx629480) diluted at 1:2000.
 Blue: Hoechst 33342 staining.
 Scale bar = 10 μm.


GTx629480 WB Image

Factor X antibody [GT268] detects Factor X protein by Western blot analysis.
 A. 30 μg 293T whole cell lysate/extract
 B. 30 μg whole cell lysate/extract of human F10-transfected 293T cells
 10 % SDS-PAGE
 Factor X antibody [GT268] (GTx629480) dilution: 1:10000



For full product information, images and publications, please visit our [website](https://www.genetex.com).