

PUF60 antibody [GT917]

Cat. No. GTX629887

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	1:100-1:1000

Not tested in other applications.

Calculated MW 60 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human PUF60. The exact sequence is proprietary.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX629887 IHC-P Image

PUF60 antibody [GT917] detects PUF60 protein at nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded mouse testis.

PUF60 stained by PUF60 antibody [GT917] (GTX629887) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




GTX629887 IHC-P Image

PUF60 antibody [GT917] detects PUF60 protein at nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded rat colon.

PUF60 stained by PUF60 antibody [GT917] (GTX629887) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




GTX629887 WB Image

PUF60 antibody [GT917] detects PUF60 protein by Western blot analysis.

A. 50 µg rat brain lysate/extract

7.5 % SDS-PAGE

PUF60 antibody [GT917] (GTX629887) dilution: 1:1000



GTX629887 IHC-P Image

PUF60 antibody [GT917] detects PUF60 protein at nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.

PUF60 stained by PUF60 antibody [GT917] (GTX629887) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).