

p39 antibody [EPR5074]

Cat. No. GTX62997

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P, IP
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IP	Assay dependent

Not tested in other applications.


Calculated MW 39 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Buffer, 0.05% BSA, 50% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) corresponding to p39.
Purification	Tissue culture supernatant
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. RabMAb® technology is covered by the following U.S. Patents, No. 5,675,063 and/or 7,429,488.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx62997 WB Image

WB analysis of lysates from (1) human fetal brain, (2) SH-SY5Y, (3) mouse brain and (4) rat brain (10 µg per lane) using p39 antibody [EPR5074] at a dilution of 1:1,000.



GTx62997 IHC-P Image

IHC-P analysis of human brain tissue using p39 antibody [EPR5074] at a dilution of 1:100.



For full product information, images and publications, please visit our [website](#).