

LIN7A antibody [GT879]

Cat. No. GTX631107

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IP
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	1:100-1:1000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 26 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human LIN7A. The exact sequence is proprietary.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX631107 IP Image

Immunoprecipitation of LIN7A protein from MDA-MB-231 whole cell extracts using 5 µg of LIN7A antibody [GT879] (GTX631107).

Western blot analysis was performed using LIN7A antibody [GT879] (GTX631107).

EasyBlot anti-Rabbit IgG (GTX221666-01) was used as a secondary reagent.



GTX631107 ICC/IF Image

LIN7A antibody [GT879] detects LIN7A protein at junction by immunofluorescent analysis.

Sample: MCF-7 cells were fixed in ice-cold MeOH for 5 min.

Green: LIN7A protein stained by LIN7A antibody [GT879] (GTX631107) diluted at 1:500.

Blue: Hoechst 33342 staining.


GTX631107 WB Image

LIN7A antibody [GT879] detects LIN7A protein by western blot analysis.

A. 30 µg MCF-7 whole cell extract

B. 30 µgMDA-MB-231 whole cell extract

12 % SDS-PAGE

LIN7A antibody [GT879] (GTX631107) dilution: 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).