

58K Golgi protein antibody [GT575]

Cat. No. GTX631588

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, IP
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IP	1:100-1:500

Not tested in other applications.

Calculated MW 59 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human 58K Golgi protein. The exact sequence is proprietary.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX631588 IP Image


Immunoprecipitation of FTCD protein from HepG2 whole cell extracts using 5 µg of FTCD antibody [GT575] (GTX631588).

Western blot analysis was performed using FTCD antibody [GT575] (GTX631588).

EasyBlot anti-Mouse IgG (GTX221667-01) was used as a secondary reagent.


GTX631588 WB Image

FTCD antibody detects FTCD protein by western blot analysis. Mouse tissue extracts (50 µg) was separated by 7.5 % SDS-PAGE, and the membrane was blotted with FTCD antibody (GTX631588) diluted by 1:500.


GTX631588 WB Image

FTCD antibody detects FTCD protein by western blot analysis. Whole cell extracts (30 µg) was separated by 7.5 % SDS-PAGE, and blotted with FTCD antibody (GTX631588) diluted by 1:1000.



For full product information, images and publications, please visit our [website](https://www.genetex.com).