

ARID5A antibody [GT921]

Cat. No. GTX631938

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, FCM
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:10000
FCM	1:20-1:200

Not tested in other applications.

Observed MW (kDa) 72 kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein encompassing a sequence within the center region of human ARID5A. The exact sequence is proprietary.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX631938 WB Image

THP-1 whole cell and nuclear extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with ARID5A antibody [GT921] (GTX631938) diluted at 1:500. The signal was developed with Trident ECL plus-Enhanced.


GTX631938 FCM Image


ARID5A antibody [GT921] (GTX631938) detects ARID5A protein by flow cytometry analysis.

Sample: THP-1 cell.

Black: Unlabelled sample was used as a control.

Red: ARID5A antibody [GT921] (GTX631938) dilution: 1:20.

Acquisition of 20,000 events were collected for FACS analysis.


GTX631938 WB Image

ARID5A antibody detects ARID5A protein by western blot analysis. Non-transfected (-) and ARID5A-transfected (+, including 3xFlag-tag) 293T whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with ARID5A antibody (GTX631938) diluted by 1:5000.



For full product information, images and publications, please visit our [website](https://www.genetex.com).