

TOM1L1 antibody [3F12]

Cat. No. GTX632246

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2.34 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen used to generate this antibody corresponds to human TOM1L1.
Purification	Purified by antigen-affinity chromatography.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx632246 ICC/IF Image

ICC/IF analysis of EGF stimulated A431 cells using GTx632246 TOM1L1 antibody [3F12].

Red : Primary antibody

Green : Polyclonal TOM1L1 antibody


Yellow : Merged color of red and green


GTx632246 WB Image

WB analysis of HeLa lysate transfected with TOM1L1 using GTx632246 TOM1L1 antibody [3F12].

T : Transfected HeLa with HA-TOM1L1 fL

L : HeLa lysate


GTx632246 WB Image

TOM1L1 antibody detects TOM1L1 protein by western blot analysis. Non-transfected (-) and TOM1L1-transfected (+) 293T whole cell extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with TOM1L1 antibody (GTx632246) diluted at 1:1000.



For full product information, images and publications, please visit our [website](https://www.genetex.com).