

NUF2 antibody [GT312]

Cat. No. GTX632247

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-P, IP
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	Assay dependent
IP	1:100-1:500

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein corresponding to human NUF2
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX632247 IHC-P Image

NUF2 antibody [GT312] detects NUF2 protein at nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded human breast carcinoma.


NUF2 stained by NUF2 antibody [GT312] (GTX632247) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX632247 WB Image

HeLa whole cell and nuclear extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with NUF2 antibody [GT312] (GTX632247) diluted at 1:1000.



GTX632247 IP Image

Immunoprecipitation of NUF2 protein from HeLa whole cell extracts using 5 µg of NUF2 antibody [GT312] (GTX632247) .

Western blot analysis was performed using NUF2 antibody [GT312] (GTX632247) diluted at 1:500.

EasyBlot anti-Mouse IgG (GTX221667-01) was used as a secondary reagent.



For full product information, images and publications, please visit our [website](https://www.genetex.com).