

AGTR1 antibody [EPR3873]

Cat. No. GTX63229

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µl

Applications

Application Note

Recommended Starting Dilutions:

For WB: Use at 1:1,000 - 10,000

Not yet tested in other applications. Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 41 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	50mM Tris-Glycine, 150mM NaCl, 0.05% BSA, 40% Glycerol
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	A synthetic peptide corresponding to residues on the N-terminus in human Type-1 Angiotensin II receptor was used as an immunogen.
Purification	Tissue culture supernatant
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. RabMAb® technology is covered by the following U.S. Patents, No. 5,675,063 and/or 7,429,488.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX63229 WB Image

Western blot analysis on (A) C2C12, (B) PC-12, (C) HeLa, and (D) HepG2 cell lysates using anti-Type-1 Angiotensin II receptor RabMAb (cat. #GTX63229).



For full product information, images and publications, please visit our [website](https://www.genetex.com).