

## ASFV p54 antibody [HL1218]

**Cat. No. GTX636562**

<b>Host</b>	Rabbit
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB, ICC/IF, ELISA, Lateral Flow, Sandwich ELISA
<b>Reactivity</b>	African swine fever virus

**Package**

100 µl, 25 µl

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:5000-1:500000
ICC/IF	Assay dependent
ELISA	Assay dependent
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture : GTX635690 / GTX635691 / GTX636563 / GTX636700 / GTX636701, Detection : GTX636562****Capture : GTX636562, Detection : GTX636563.****Please notice that GTX636563 needs to be conjugated to HRP when paired with GTX636562. Please contact us for custom HRP-conjugated antibody.**

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant protein encompassing a sequence within the center region of ASFV p54. (African swine fever virus (strain: BA71V))
<b>Purification</b>	Affinity purified by Protein A.
<b>Conjugation</b>	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

#### Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

#### DATA IMAGES




#### GTx636562 ICC/IF Image

ASFV p54 antibody [HL1218] detects ASFV p54 protein by immunofluorescent analysis.  
Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.  
Green: ASFV p54 stained by ASFV p54 antibody [HL1218] (GTx636562) diluted at 1:500.  
Blue: Fluoroshield with DAPI (GTx30920).



#### GTx636562 WB Image

Non-transfected (–) and transfected (+) 293T whole cell extracts (30 μg) were separated by 15% SDS-PAGE, and the membrane was blotted with ASFV p54 antibody [HL1218] (GTx636562) diluted at 1:100000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.




#### GTx636562 ELISA Image

Sandwich ELISA detection of recombinant ASFV p54 using antibodies as below.

**Capture:** ASFV p54 antibody [HL1218] (GTx636562) (5 μg/mL)

**Detection:** HRP-conjugated ASFV p54 antibody [HL1219] (GTx636563) (0.1 μg/mL)

Please notice that GTx636563 needs to be conjugated to HRP to function as the detection antibody when paired with GTx636562. Please contact us for custom HRP-conjugated antibody.



#### GTx636562 Lateral Flow Image

Detection of HEK293 expressed, ASFV p54 (ECD) protein by lateral flow assay using the monoclonal antibody pair.

**Capture:** ASFV p54 antibody (GTx635690 [GT853])

**Detection:** ASFV p54 antibody (GTx636562 [HL1218])

**Samples (40 ng) :**

1. ASFV p54 (ECD) protein (GTx137950-pro)
2. Lysis buffer



For full product information, images and publications, please visit our [website](https://www.genetex.com).