

EPO antibody [HL1794] (HRP)

Cat. No. GTX637446-01

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ELISA
Reactivity	Human

Package

50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Note : Capture : GTX637445, Detection : GTX637446-01

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.83 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human EPO.
Purification	Affinity purified by Protein A.
Conjugation	Horseradish peroxidase(HRP)

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES



GTX637446-01 ELISA Image

Sandwich ELISA detection of recombinant (CHO cells expressed, full-length human EPO protein) using antibodies as below.

Capture: EPO antibody [HL1793] (GTX637445) (5 µg/mL)

Detection: HRP-conjugated EPO antibody [HL1794] (GTX637446-01) (1 µg/mL)



For full product information, images and publications, please visit our [website](#).