

FLI1 antibody [HL1805]

Cat. No. GTX637457

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:5000-1:100000
IP	Assay dependent

Not tested in other applications.

Observed MW (kDa) 51 kDa.

Product Note This antibody is specific for human FLI1 protein, and it does not cross react with human ERG protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human FLI1
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX637457 WB Image

Non-transfected (–) and transfected (+) 293T whole cell extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with FLI1 antibody [HL1805] (GTX637457) diluted at 1:100000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.


GTX637457 WB Image

Various whole cell extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with FLI1 antibody [HL1805] (GTX637457) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.


GTX637457 IP Image

Immunoprecipitation of FLI1 protein from Raji whole cell extract using 5 µg of FLI1 antibody [HL1805] (GTX637457).

Western blot analysis was performed using FLI1 antibody [HL1805] (GTX637457).

EasyBlot HRP-conjugated anti rabbit IgG antibody (GTX221666-01).



For full product information, images and publications, please visit our [website](https://www.genetex.com).