

Eosinophil derived neurotoxin antibody [HL2166]

Cat. No. GTX638144

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Observed MW (kDa) 17 (de-glycosylated form), 18-40 (glycosylated form) kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human Eosinophil derived neurotoxin
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX638144 WB Image

Untreated (–) and treated (+) THP-1 whole cell extracts (30 µg) were separated by 15% SDS-PAGE, and the membrane was blotted with Eosinophil derived neurotoxin antibody [HL2166] (GTX638144) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX638144 WB Image

Untreated (–) and treated (+) THP-1 whole cell extract (30 µg) were separated by 15% SDS-PAGE, and the membrane was blotted with Eosinophil derived neurotoxin antibody [HL2166] (GTX638144) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX638144 IHC-P Image

Eosinophil derived neurotoxin antibody [HL2166] detects Eosinophil derived neurotoxin protein at cytoplasm by immunohistochemical analysis.

Sample: Paraffin-embedded human kidney.

Eosinophil derived neurotoxin stained by Eosinophil derived neurotoxin antibody [HL2166] (GTX638144) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).