

Chikungunya virus nsP4 antibody [HL2580]

Cat. No. GTX638958

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Chikungunya Virus

References (2)

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1:100-1:1000

Not tested in other applications.

Properties


Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of Chikungunya virus nsP4 (strain S27-African prototype)
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX638958 ICC/IF Image

Immunofluorescent analysis of Chikungunya virus infected cells using Chikungunya virus nsP4 antibody [HL2580] antibody (GTX638958).

Sample: Chikungunya virus infected cells slide.

Green: Chikungunya virus nsP4 antibody [HL2580] antibody (GTX638958) diluted at 1:100.

Blue: Fluoroshield with DAPI (GTX30920).



GTX638958 ICC/IF Image

Chikungunya virus nsP4 antibody [HL2580] detects Chikungunya virus nsP4 protein by immunofluorescent analysis.


Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: Chikungunya virus nsP4 stained by Chikungunya virus nsP4 antibody [HL2580] (GTX638958) diluted at 1:500.

Blue: Fluoroshield with DAPI (GTX30920).


GTX638958 WB Image

Non-transfected (–) and transfected (+) 293T whole cell extracts (30 μg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with Chikungunya virus nsP4 antibody [HL2580] (GTX638958) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.


GTX638958 WB Image

Chikungunya viral lysate (0.1 μg) was separated by 7.5% SDS-PAGE, and the membrane was blotted with Chikungunya virus nsP4 antibody [HL2580] (GTX638958) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.



For full product information, images and publications, please visit our [website](https://www.genetex.com).