

Fibrillarin + FBLL1 antibody [HL2588]

Cat. No. GTX639039

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package

100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1:100-1:1000
IHC-P	Assay dependent

Not tested in other applications.

Observed MW (kDa) 35-38 kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human fibrillarin.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx639039 ICC/IF Image


Fibrillarin + FBLL1 antibody [HL2588] detects Fibrillarin + FBLL1 protein at nuclear envelope and nucleus by immunofluorescent analysis.

Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: Fibrillarin + FBLL1 stained by Fibrillarin + FBLL1 antibody [HL2588] (GTx639039) diluted at 1:500.


Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (GTx628802) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTx30920).




GTx639039 WB Image

Various whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with Fibrillarin + FBLL1 antibody [HL2588] (GTx639039) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.



GTx639039 WB Image

Various whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with Fibrillarin + FBLL1 antibody [HL2588] (GTx639039) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.




GTx639039 IHC-P Image

Fibrillarin + FBLL1 antibody [HL2588] detects Fibrillarin + FBLL1 protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat tissues.

Fibrillarin + FBLL1 stained by Fibrillarin + FBLL1 antibody [HL2588] (GTx639039) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTx639039 IHC-P Image

Fibrillarin + FBLL1 antibody [HL2588] detects Fibrillarin + FBLL1 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

Fibrillarin + FBLL1 stained by Fibrillarin + FBLL1 antibody [HL2588] (GTx639039) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).