

GPX8 antibody [HL2635]

Cat. No. GTX639091

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat, Dog

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	Assay dependent

Not tested in other applications.

Observed MW (kDa) 24 kDa.

Product Note This antibody may not be suitable for mouse / rat WB assay, as its mouse / rat reactivity is based on IHC-P testing.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human GPX8.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX639091 IHC-P Image

GPX8 antibody [HL2635] detects GPX8 protein by immunohistochemical analysis.

Sample: Paraffin-embedded rat tissues.

GPX8 stained by GPX8 antibody [HL2635] (GTX639091) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




GTX639091 IHC-P Image

GPX8 antibody [HL2635] detects GPX8 protein at cell membrane by immunohistochemical analysis.

Sample: Paraffin-embedded mouse kidney.


GPX8 stained by GPX8 antibody [HL2635] (GTX639091) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX639091 WB Image

Various whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with GPX8 antibody [HL2635] (GTX639091) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.



GTX639091 WB Image

Whole cell extract (30 µg) was separated by 12% SDS-PAGE, and the membrane was blotted with GPX8 antibody [HL2635] (GTX639091) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).