

Glucocorticoid Receptor antibody [HMV304] HistoMAX™

Cat. No. GTX639927

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

References (1)

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)

[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution Recommended dilution

IHC-P 1:100-1:200

Note : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in Antigen retrieval buffer, pH 7.8 (GTX443621).

Not tested in other applications.

Product Note **Highly recommended for IHC-P in human tissues.**

Properties

Form Liquid

Buffer PBS, 2% BSA

Preservative 0.025% ProClin 300

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 56 µg/ml (Please refer to the vial label for the specific concentration.)

Immunogen Recombinant protein encompassing a sequence within the center region of human Glucocorticoid Receptor. The exact sequence is proprietary.

Purification Affinity purified by Protein A.

Conjugation Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 13 Page 1 of 2


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX639927 IHC-P Image

IHC-P analysis of human fallopian tube mucosa tissue using GTX639927 Glucocorticoid Receptor antibody [HMV304] HistoMAX™.
Strong glucocorticoid receptor positivity of all epithelial cells.


GTX639927 IHC-P Image

IHC-P analysis of human adrenal gland tissue using GTX639927 Glucocorticoid Receptor antibody [HMV304] HistoMAX™.
Moderate to strong GR positivity of medullary cells while adrenocortical cells are largely GR negative.


GTX639927 IHC-P Image

IHC-P analysis of human high-grade serous ovarian cancer (HGSOC) tissue using GTX639927 Glucocorticoid Receptor antibody [HMV304] HistoMAX™.
A strong glucocorticoid receptor immunostaining of all tumor cells.



For full product information, images and publications, please visit our [website](https://www.genetex.com).