

MEF2C antibody [HL2811]

Cat. No. GTX639952

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Observed MW (kDa) 43-60 kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human MEF2C.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX639952 WB Image

Various whole cell extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with MEF2C antibody [HL2811] (GTX639952) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.




GTX639952 IHC-P Image

MEF2C antibody [HL2811] detects MEF2C protein at nucleus by immunohistochemical analysis.
Sample: Paraffin-embedded rat heart.

MEF2C stained by MEF2C antibody [HL2811] (GTX639952) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




GTX639952 IHC-P Image

MEF2C antibody [HL2811] detects MEF2C protein at nucleus by immunohistochemical analysis.

Sample: Paraffin-embedded mouse brain.

MEF2C stained by MEF2C antibody [HL2811] (GTX639952) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min




GTX639952 ICC/IF Image

MEF2C antibody [HL2811] detects MEF2C protein at nucleus by immunofluorescent analysis.

Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: MEF2C stained by MEF2C antibody [HL2811] (GTX639952) diluted at 1:500.

Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (GTX628802) diluted at 1:1000.



GTX639952 WB Image

Various whole cell extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with MEF2C antibody [HL2811] (GTX639952) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.



For full product information, images and publications, please visit our [website](#).