

C3a Receptor antibody [HL2855]

Cat. No. GTX640126

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF, IHC-P
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100-1:1000
IHC-P	Assay dependent

Not tested in other applications.

Product Note This antibody was raised against human C3a Receptor Extracellular domain.

Properties


Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human C3a Receptor. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX640126 ICC/IF Image

C3a Receptor antibody [HL2855] detects C3a Receptor protein by immunofluorescent analysis.

Sample: Mock and transfected 293T cells were fixed in ice-cold MeOH for 5 min.

Green: C3a Receptor stained by C3a Receptor antibody [HL2855] (GTX640126) diluted at 1:500.

Blue: Fluoroshield with DAPI (GTX30920).


GTX640126 IHC-P Image

C3a Receptor antibody [HL2855] detects C3a Receptor protein at cell membrane by immunohistochemical analysis.

Sample: Paraffin-embedded human colon.

C3a Receptor stained by C3a Receptor antibody [HL2855] (GTX640126) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).