

mGluR5 antibody [HL2876]

Cat. No. GTX640153

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent

Not tested in other applications.

Observed MW (kDa) 130 (monomer), > 250 (dimer or multimer) kDa.

Product Note This antibody was raised against human mGluR5 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain human of mGluR5. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX640153 WB Image


Various tissue extracts (50 μ g) were separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL2876] (GTX640153) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.


GTX640153 WB Image


Rat tissue extract (50 μ g) was separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL2877] (GTX640153) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.


GTX640153 WB Image

Boiled and unboiled various tissue extracts (50 μ g) were separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL2876] (GTX640153) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.


GTX640153 WB Image

Boiled and unboiled human tissue extract (50 μ g) were separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL2876] (GTX640153) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident femto Western HRP Substrate.


GTX640153 ICC/IF Image

mGluR5 antibody [HL2876] detects mGluR5 protein by immunofluorescent analysis. Sample: DIV11 rat E18 hippocampal neurons and glial cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: mGluR5 stained by mGluR5 antibody [HL2876] (GTX640153) diluted at 1:250. Red: Tau, an axon marker, stained by Tau antibody [GT287] (GTX634809) diluted at 1:500. Blue: Fluoroshield with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](https://www.genetex.com).