

beta 2 Adrenergic Receptor antibody [HL3099]

Cat. No. GTX640545

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	FCM
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	1:50-1:200

Note : Do not need a permeabilization step.

Not tested in other applications.

Product Note This antibody was raised against human beta 2 Adrenergic Receptor Extracellular domain.


Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human beta 2 Adrenergic Receptor. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note
For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx640545 FCM Image

beta 2 Adrenergic Receptor antibody [HL3099] (GTx640545) detects ADRB2 protein by flow cytometry analysis.

Sample: non-fixed and non-permeabilized ACHN or HeLa cells were stained with 5 µg/ml GTx640545 (red) or a Rabbit monoclonal IgG isotype control (gray) in 1% FBS/ 5 mM EDTA/PBS at 4°C for 1 hour.



For full product information, images and publications, please visit our [website](https://www.genetex.com).