

CA19-9 antibody [HMV333] HistoMAX™

Cat. No. GTX640565

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Species independent

Package

500 µl, 100 µl

PRODUCT

This antibody was validated on 76 different Normal Tissues by IHC-P.

Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

IHC-P	1:100-1:200
-------	-------------

Note : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in Antigen retrieval buffer, pH 7.8 (GTX443621).

Not tested in other applications.

Product Note

Highly recommended for IHC-P in human tissues.

Properties

Form	Liquid
Buffer	PBS, 2% BSA, 20% Glycerol
Preservative	0.025% ProClin 300
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.22 mg/ml
Immunogen	The immunogen used to generate this antibody corresponds to CA19-9.
Purification	Affinity purified by Protein G.
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 09 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX640565 IHC-P Image

IHC-P analysis of human esophagus tissue using GTX640565 CA19-9 antibody [HMOV333] HistoMAX™. Strong membranous CA19-9 staining of a significant fraction of cells in superficial layers of the squamous epithelium.



GTX640565 IHC-P Image

IHC-P analysis of human mucosa from appendix tissue using GTX640565 CA19-9 antibody [HMOV333] HistoMAX™. Moderate to strong apical membranous CA19-9 staining of most superficial epithelial cells. In addition, strong membranous and cytoplasmic positivity of a large fraction of goblet cells.



GTX640565 IHC-P Image

IHC-P analysis of human pancreas tissue using GTX640565 CA19-9 antibody [HMOV333] HistoMAX™. A CA19-9 negative acinar cell carcinoma.



For full product information, images and publications, please visit our [website](#).