

FSH Receptor antibody [HL3233]

Cat. No. GTX640874

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, FCM
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
FCM	Assay dependent

Not tested in other applications.

Observed MW (kDa) 70 kDa.

Product Note This antibody was raised against human FSH Receptor Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human FSH Receptor. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX640874 WB Image

Human tissue extract (10 µg) was separated by 7.5% SDS-PAGE, and the membrane was blotted with FSH Receptor antibody [HL3233] (GTX640874) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX640874 FCM Image

FSH Receptor antibody [HL3233] (GTX640874) detects FSH Receptor protein by flow cytometry analysis. Sample: Non-transfected and transfected 293T cells were fixed in 4% paraformaldehyde at 4°C for 15 min and washed with 0.1% saponin buffer. The cells were following stained with FSH Receptor antibody [HL3233] (GTX640874) diluted at 1:50 or a Rabbit monoclonal IgG (GTX35035) isotype control at 4°C for 1 hour.



For full product information, images and publications, please visit our [website](https://www.genetex.com).