

Cnd2 (S. pombe) antibody

Cat. No. GTX64102

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Schizosaccharomyces pombe

Package
100 µl

Applications

Application Note

Western blotting (500 fold dilution)

Calculated MW	83 kDa. (Note)
Product Note	Reacts with S. pombe Cnd2 protein. Not tested for other species

Properties

Form	Liquid
Buffer	Serum
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant GST-Cnd2 fusion protein expressed in E. coli
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES


GTX64102 WB Image

Identification of Cnd2 protein and Cnd2-GFP protein in crude cell extract by anti-Cnd2 antibody at 1/500 dilution.


GTX64102 WB Image

Identification of Cnd2 protein in the crude extract of *S. pombe* by western blotting with this antiserum.


GTX64102 ICC/IF Image

Detection of Cnd2 protein by immunofluorescent staining of *S. pombe* cell.



For full product information, images and publications, please visit our [website](https://www.genetex.com).