

Chromogranin C antibody [HL3311]

Cat. No. GTX641073

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	Assay dependent
IHC-Fr	Assay dependent

Not tested in other applications.

Observed MW (kDa)	88-95 kDa.
Product Note	Highly recommended for IHC-P in human and mouse tissues.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human Chromogranin C.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX641073 WB Image

SH-SY5Y whole cell extract and conditioned medium (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with Chromogranin C antibody [HL3311] (GTX641073) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX641073 WB Image

Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with Chromogranin C antibody [HL3311] (GTX641073) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.




GTX641073 IHC-P Image

Chromogranin C antibody [HL3311] detects Chromogranin C protein by immunohistochemical analysis. Sample: Paraffin-embedded rat eye.

Chromogranin C stained by Chromogranin C antibody [HL3311] (GTX641073) diluted at 1:250.

Red: beta Tubulin 3/ Tuj1, stained by beta Tubulin 3/ Tuj1 antibody [GT11710] (GTX631836) diluted at 1:500.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



GTX641073 IHC-P Image

Chromogranin C antibody [HL3311] detects Chromogranin C protein by immunohistochemical analysis. Sample: Paraffin-embedded human tissues.

Chromogranin C stained by Chromogranin C antibody [HL3311] (GTX641073) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min

Corresponding RNA expression data for the same tissues are based on NCBI.



For full product information, images and publications, please visit our [website](https://www.genetex.com).