

Suppressor of Fused antibody [HL3374]

Cat. No. GTX641194

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Observed MW (kDa) 54 kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml
Immunogen	Recombinant fragment of human Suppressor of Fused.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX641194 WB Image

Rat whole cell extracts (30 µg) was separated by 10% SDS-PAGE, and the membrane was blotted with Suppressor of Fused antibody [HL3374] (GTX641194) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX641194 WB Image

Various whole cell extracts (30 µg) were separated by 10% SDS-PAGE, and the membrane was blotted with Suppressor of Fused antibody [HL3374] (GTX641194) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



GTX641194 IHC-P Image


Suppressor of Fused antibody [HL3374] detects Suppressor of Fused protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse tissues.

Suppressor of Fused stained by Suppressor of Fused antibody [HL3374] (GTX641194) diluted at 1:100.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min

Corresponding RNA expression data for the same tissues are based on NCBI.



GTX641194 IHC-P Image

Suppressor of Fused antibody [HL3374] detects Suppressor of Fused protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tissues.

Suppressor of Fused stained by Suppressor of Fused antibody [HL3374] (GTX641194) diluted at 1:250.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min

Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.



For full product information, images and publications, please visit our [website](https://www.genetex.com).