

HKDC1 antibody [HL3464]

Cat. No. GTX641349

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
ICC/IF	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Observed MW (kDa) 120 kDa.

Product Note Highly recommended for IHC-P in mouse tissues.
This antibody is specific for human HKDC1 protein, and it does not cross react with human HK1, HK2, and HK3 protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human HKDC1.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx641349 IHC-P Image


HKDC1 antibody [HL3464] detects HKDC1 protein by immunohistochemical analysis.

Sample: Paraffin-embedded human tissues.


HKDC1 stained by HKDC1 antibody [HL3464] (GTx641349) diluted at 1:200.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.


GTx641349 WB Image

Various tissue extracts (50 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with HKDC1 antibody [HL3464] (GTx641349) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx641349 WB Image

Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with HKDC1 antibody [HL3464] (GTx641349) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.


GTx641349 ICC/IF Image

HKDC1 antibody [HL3464] detects HKDC1 protein at mitochondria by immunofluorescent analysis.

Sample: AsPC-1 cells were fixed in ice-cold MeOH for 5 min.

Green: HKDC1 stained by HKDC1 antibody [HL3464] (GTx641349) diluted at 1:500.

Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (GTx628802) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTx30920).



For full product information, images and publications, please visit our [website](https://www.genetex.com).