

beta Catenin (phospho Ser675) antibody [HL3471]

Cat. No. GTX641356

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
Not tested in other applications.	

Observed MW (kDa) 92 kDa.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Carrier-protein conjugated synthetic peptide surrounding phospho Ser65 of human beta Catenin. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX641356 WB Image

Untreated (-) and treated (+) HeLa whole cell extract (30 µg) were separated by 7.5% SDS-PAGE, and the membranes were blotted with beta Catenin (phospho Ser675) antibody [HL3471] (GTX641356) diluted at 1:1000 and beta Catenin antibody (GTX637403) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident femto Western HRP Substrate.


GTX641356 WB Image

Untreated (-) and treated (+) NIH-3T3 whole cell extract (30 µg) were separated by 7.5% SDS-PAGE, and the membranes were blotted with beta Catenin (phospho Ser675) antibody [HL3471] (GTX641356) diluted at 1:1000 and beta Catenin antibody (GTX637403) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).