

Mycoplasma pneumoniae P1 protein antibody [GT50]

Cat. No. GTX641374

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	ELISA, Sandwich ELISA
Reactivity	Mycoplasma pneumoniae

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX641376, Detection : GTX640845

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human Mycoplasma pneumoniae P1 protein.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX641374 ELISA Image

Sandwich ELISA detection of recombinant E. coli expressed, partial Mycoplasma pneumoniae P1 protein using antibodies as below.

Capture: Mycoplasma pneumoniae P1 protein antibody [GT50] (GTX641374) (5 µg/mL)


Detection: Mycoplasma pneumoniae P1 protein anti body [HL3204] (GTX640844) (1 µg/mL)


GTX641374 ELISA Image

Sandwich ELISA detection of recombinant E. coli expressed, partial Mycoplasma pneumoniae P1 protein using antibodies as below.

Capture: Mycoplasma pneumoniae P1 protein antibody [GT50] (GTX641374) (5 µg/mL)

Detection: Mycoplasma pneumoniae P1 protein anti body [HL3205] (GTX640845) (1 µg/mL)


GTX641374 ELISA Image

Sandwich ELISA detection of recombinant E. coli expressed, partial Mycoplasma pneumoniae P1 protein using antibodies as below.

Capture: Mycoplasma pneumoniae P1 protein antibody [GT50] (GTX641374) (5 µg/mL)

Detection: Mycoplasma pneumoniae P1 protein anti body [HL3203] (GTX640843) (1 µg/mL)



For full product information, images and publications, please visit our [website](https://www.genetex.com).