

mGluR5 antibody [HL3509]

Cat. No. GTX641397

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr, FCM
Reactivity	Human, Mouse, Rat

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Observed MW (kDa) 130 (monomer), > 250 (dimer or multimer) kDa.**Product Note** This antibody was raised against mGluR5 Intracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment encompassing a sequence within the Intracellular domain of human mGluR5. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX641397 IHC-P Image


mGluR5 antibody [HL3509] detects mGluR5 protein by immunohistochemical analysis.

Sample: Paraffin-embedded mouse hippocampus.


mGluR5 stained by mGluR5 antibody [HL3509] (GTX641397) diluted at 1:400.

Antigen Retrieval: Citrate buffer, pH 6.0, 15 min


*The competitor is not affiliated with GeneTex and does not endorse this product.


GTX641397 WB Image

Boiled and unboiled human tissue extract (5 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL3509] (GTX641397) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.


GTX641397 WB Image

Unboiled various tissue extracts (50 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL3509] (GTX641397) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.


GTX641397 WB Image

Unboiled various tissue extracts (50 µg) were separated by 5% SDS-PAGE, and the membrane was blotted with mGluR5 antibody [HL3509] (GTX641397) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTX213110-01) was used to detect the primary antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).