

Streptococcus NADase Inhibitor antibody

Cat. No. GTX64170

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, ELISA
Reactivity	Staphylococcus

Package
100 µl

Applications

Application Note

Western blotting (1:3,000 dilution)


Properties

Form	Liquid
Buffer	Serum
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified recombinant His-tagged NADase inhibitor of <i>S. dysgalactiae</i> expressed in <i>E. coli</i>
Purification	Unpurified
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX64170 WB Image

Identification of NADase inhibitor in crude extracts by western blotting with anti NADase inhibitor antibody. The antibody was used at 1/3,000 dilution. 1. Cell lysate of *S. dysgalactiae* 2. Cell lysate of *E. coli* overexpressing His-tagged NADase NADase inhibitor. 3. Culture supernatant of *S. dysgalactiae* 4. Molecular size markers.



For full product information, images and publications, please visit our [website](https://www.genetex.com).