

Sup35 / PSI+ (S. cerevisiae) antibody

Cat. No. GTX64173

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Saccharomyces cerevisiae

Package
100 µl

Applications

Application Note

Western blotting (1,000-2,000 fold dilution)

Calculated MW 77 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to a.a. 494-507 of Sup35
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX64173 WB Image

Identification of Sup35 protein in crude extracts of *S. cerevisiae* by western blotting with anti-Sup35 antibody 1. Endogenous expression in *S. cerevisiae*. 2. Overexpression of Sup35 protein by introduction of a plasmid carrying Sup35 gene into *S. cerevisiae*. Anti-Sup35 antibody was used at 1/2,000 dilution. Sup35 protein was identified about 77 kDa position (Molecular mass of Sup35 is 76.6 kDa)



For full product information, images and publications, please visit our [website](https://www.genetex.com).