

# Taf11 (S. cerevisiae) antibody

**Cat. No. GTX64176**

|                     |                          |
|---------------------|--------------------------|
| <b>Host</b>         | Rabbit                   |
| <b>Clonality</b>    | Polyclonal               |
| <b>Isotype</b>      | IgG                      |
| <b>Applications</b> | WB                       |
| <b>Reactivity</b>   | Saccharomyces cerevisiae |

**Package**  
100 µl

## Applications

### Application Note

Western blotting (1:500-1:1,000 dilution).

**Calculated MW** 41 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | Serum                                                                                                                                                                                                                      |
| <b>Preservative</b>  | 0.1% Sodium azide                                                                                                                                                                                                          |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| <b>Immunogen</b>     | Recombinant His-tagged S. cerevisiae Taf11 (1-176 aa) expressed in E. coli.                                                                                                                                                |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX64176 WB Image

Detection of endogenous Taf11 in whole cell extract of *S. cerevisiae* by Western blotting, using the Taf11 antibody. The antiserum was used at 1/500 dilution.



For full product information, images and publications, please visit our [website](https://www.genetex.com).