

Taf4 (S. cerevisiae) antibody

Cat. No. GTX64179

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Saccharomyces cerevisiae

Package
100 µl

Applications

Application Note

Western blotting. (1:1,000-1:5,000 dilution)

Calculated MW 42 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant His-tagged N-terminal domain of Taf4p protein (1-200aa)
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx64179 WB Image

Detection of Taf4p by Western blotting Lane 1: Extract of budding yeast The antiserum was diluted 5,000 fold before use. The molecular mass of Taf4 is 48 kDa



For full product information, images and publications, please visit our [website](https://www.genetex.com).