

Taf6 (S. cerevisiae) antibody

Cat. No. GTX64180

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Saccharomyces cerevisiae

Package
100 µl

Applications

Application Note

Western blotting (1:1,000)

Calculated MW 58 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant His tagged Taf6 protein (1-200 aa)
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX64180 WB Image

Identification of endogenous Taf6 protein in crude cell extract of *S. cerevisiae*. The anti-Taf6 antibody was used at 1/1,000 dilution. Molecular mass of Taf6 is 60 kDa



For full product information, images and publications, please visit our [website](https://www.genetex.com).