

Toa1 (S. cerevisiae) antibody

Cat. No. GTX64185

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Saccharomyces cerevisiae

Package
100 µl

Applications

Application Note

Western blotting (1:1,000-1:5,000)

Calculated MW 32 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Recombinant His-tagged Toa1 protein (1-286aa) produced in E. coli
Purification	Unpurified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX64185 WB Image

Detection of endogenous Toa1p by Western blotting Lane M: Protein size markers (kDa) Lane1: Extract of budding yeast The antiserum was used at 1/5,000 dilution Molecular mass of Toa1 is 32.2 kDa, but it shows apparent size of 45~47 kDa in SDS-PAGE-WB as in other publications.



For full product information, images and publications, please visit our [website](https://www.genetex.com).