

BASP1 antibody [HL3760]

Cat. No. GTX641959

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

 Review (1)Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:3000
IHC-P	Assay dependent

Note : Antigen retrieval: Citrate acid (pH6.0) is recommended.

Not tested in other applications.

Observed MW (kDa)	48 kDa.
Product Note	Highly recommended for IHC-P in mouse tissues.

Properties


Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment of human BASP1.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx641959 WB Image


Various tissue extracts (50 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with BASP1 antibody [HL3760] (GTx641959) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.


GTx641959 WB Image

Various whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with BASP1 antibody [HL3760] (GTx641959) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced. Corresponding RNA expression data for the same cell lines are based on Human Protein Atlas program.


GTx641959 IHC-P Image

BASP1 antibody [HL3760] detects BASP1 protein by immunohistochemical analysis.
Sample: Paraffin-embedded human tissues.
BASP1 stained by BASP1 antibody [HL3760] (GTx641959) diluted at 1:200.
Antigen Retrieval: Citrate buffer, pH 6.0, 15 min
Corresponding RNA levels (RPKM) in the tissues are based on NCBI database.


GTx641959 IHC-P Image

BASP1 antibody [HL3760] detects BASP1 protein by immunohistochemical analysis.
Sample: Paraffin-embedded mouse hippocampus.
BASP1 stained by BASP1 antibody [HL3760] (GTx641959) diluted at 1:250.
Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



For full product information, images and publications, please visit our [website](https://www.genetex.com).