

P2Y13 antibody [HL3775]

Cat. No. GTX641974

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	FCM
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent

Note : The FcR Blocking Reagent is recommended for use prior to the addition of the researcher's antibody to samples that contain human Fc receptor-expressing cells.

Do not need a permeabilization step.

Not tested in other applications.

Product Note This antibody was raised against human P2Y13 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human P2Y13. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES


GTX641974 FCM Image

P2Y13 antibody [HL3775] (GTX641974) detects P2RY13 protein by flow cytometry analysis.

Sample: Non-transfected and transfected 293T cells were fixed in 4% paraformaldehyde at 4°C for 15 min and washed with 0.1% saponin buffer to permeabilize the cell membrane.

The cells were following stained with P2Y13 antibody [HL3775] (GTX641974) diluted at 1:50 or a Rabbit IgG isotype control (GTX35035) control at 4°C for 1 hour.



For full product information, images and publications, please visit our [website](https://www.genetex.com).