

CXCR5 antibody [HL3777]

Cat. No. GTX641976

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100-1:1000

Not tested in other applications.

Product Note This antibody was raised against human CXCR5 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human CXCR5. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX641976 ICC/IF Image

CXCR5 antibody [HL3777] detects CXCR5 protein by immunofluorescent analysis.
 Sample: Mock and transfected 293T cells were fixed in 4% paraformaldehyde at RT for 15 min.
 Green: CXCR5 stained by CXCR5 antibody [HL3777] (GTX641976) diluted at 1:500.
 Blue: Fluoroshield with DAPI (GTX30920).



GTX641976 ICC/IF Image

CXCR5 antibody [HL3777] detects CXCR5 protein by immunofluorescent analysis.
 Sample: Daudi cells were fixed in 4% paraformaldehyde at RT for 15 min.
 Green: CXCR5 stained by CXCR5 antibody [HL3777] (GTX641976) diluted at 1:100.
 Blue: Fluoroshield with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](https://www.genetex.com).