

PAR3 antibody [HL3828]

Cat. No. GTX642052

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100-1:1000

Not tested in other applications.

Product Note This antibody was raised against human PAR3 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human PAR3. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES


GTX642052 ICC/IF Image

PAR3 antibody [HL3828] detects PAR3 protein by immunofluorescent analysis.
Sample: HEL 92.1.7 cells were fixed in 4% paraformaldehyde at RT for 15 min.
Green: PAR3 stained by PAR3 antibody [HL3828] (GTX642052) diluted at 1:500.
Blue: Fluoroshield with DAPI (GTX30920).


GTX642052 ICC/IF Image

PAR3 antibody [HL3828] detects PAR3 protein by immunofluorescent analysis.
Sample: Mock and transfected 293T cells were fixed in ice-cold MeOH for 5 min.
Green: PAR3 stained by PAR3 antibody [HL3828] (GTX642052) diluted at 1:500.
Blue: Fluoroshield with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](#).