

Adenosine A2b Receptor antibody [HL3959]

Cat. No. GTX642209

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	IHC-Fr
Reactivity	Mouse

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	1:100-1:1000

Not tested in other applications.

Product Note This antibody was raised against human Adenosine A2b Receptor Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human Adenosine A2b Receptor. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES


GTX642209 IHC-Fr Image


Adenosine A2b Receptor antibody [HL3959] detects Adenosine A2b Receptor protein by immunohistochemical analysis.

Sample: Frozen-sectioned mouse spinal cord.

Green: Adenosine A2b Receptor stained by Adenosine A2b Receptor antibody [HL3959] (GTX642209) diluted at 1:250.

Red: Vimentin antibody [GT812] (GTX629744) diluted at 1:500.

Blue: Hoechst 33342 staining.


GTX642209 IHC-Fr Image

Adenosine A2b Receptor antibody [HL3959] detects Adenosine A2b Receptor protein by immunohistochemical analysis.

Sample: Frozen-sectioned mouse duodenum.

Green: Adenosine A2b Receptor stained by Adenosine A2b Receptor antibody [HL3959] (GTX642209) diluted at 1:250.

Red: Vimentin antibody [GT812] (GTX629744) diluted at 1:500.

Blue: Hoechst 33342 staining.



For full product information, images and publications, please visit our [website](#).