

## Kiss1 antibody [HMV3972] HistoMAX™

Cat. No. GTX642298

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Monoclonal |
| Isotype      | IgG        |
| Applications | IHC-P      |
| Reactivity   | Human      |

## Package

500 µl, 100 µl

## PRODUCT

## Summary

[Go to Normal Tissue Gallery](#)[Go to Cancer Tissue Gallery](#)

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

## Suggested dilution

## Recommended dilution

IHC-P

1:100-1:200

**Note : Heat-induced antigen retrieval for 5 minutes in an autoclave at 121°C in pH 7.8 Target Retrieval Solution buffer.**

Not tested in other applications.

## Product Note

**Highly recommended for IHC-P in human tissues.**

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 2% BSA, 20% Glycerol                                                                                                                                                                                                  |
| Preservative  | 0.025% ProClin 300                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 56 µg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                  |
| Immunogen     | Recombinant fragment of human Kiss1.                                                                                                                                                                                       |
| Purification  | Affinity purified by Protein A.                                                                                                                                                                                            |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 13 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

#### DATA IMAGES




##### **GTx642298 IHC-P Image**

IHC-P analysis of human placenta (mature) tissue using GTx642298 Kiss1 antibody [HmV3972] HistoMAX™. A strong cytoplasmic KISS1 staining of the syncytiotrophoblast.



##### **GTx642298 IHC-P Image**

IHC-P analysis of first-trimester human placenta tissue using GTx642298 Kiss1 antibody [HmV3972] HistoMAX™. Intense cytoplasmic KISS1 staining of the syncytiotrophoblast while cytotrophoblast and other cells remain KISS1 negative.



##### **GTx642298 IHC-P Image**

IHC-P analysis of human placenta (mature) tissue using GTx642298 Kiss1 antibody [HmV3972] HistoMAX™. A strong cytoplasmic KISS1 staining of the majority of the syncytiotrophoblast cells.



For full product information, images and publications, please visit our [website](https://www.genetex.com).