

Adenosine A1 Receptor antibody [HL3987]

Cat. No. GTX642313

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:100-1:1000

Not tested in other applications.

Product Note This antibody was raised against human Adenosine A1 Receptor Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human Adenosine A1 Receptor. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX642313 ICC/IF Image

Adenosine A1 Receptor antibody [HL3987] detects Adenosine A1 Receptor protein by immunofluorescent analysis.

Sample: HepG2 cells were fixed in ice-cold MeOH for 5 min.

Green: Adenosine A1 Receptor stained by Adenosine A1 Receptor antibody [HL3987] (GTX642313) diluted at 1:500.

Blue: Fluoroshield with DAPI (GTX30920).



GTX642313 ICC/IF Image

Adenosine A1 Receptor antibody [HL3987] detects Adenosine A1 Receptor protein by immunofluorescent analysis.

Sample: Mock and transfected 293T cells were fixed in ice-cold MeOH for 5 min.

Green: Adenosine A1 Receptor stained by Adenosine A1 Receptor antibody [HL3987] (GTX642313) diluted at 1:500.

Blue: Fluoroshield with DAPI (GTX30920).


GTX642313 ICC/IF Image

Adenosine A1 Receptor antibody [HL3987] detects Adenosine A1 Receptor protein by immunofluorescent analysis.

Sample: ACHN cells were fixed in ice-cold MeOH for 5 min.

Green: Adenosine A1 Receptor stained by Adenosine A1 Receptor antibody [HL3987] (GTX642313) diluted at 1:500.

Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] (GTX628802) diluted at 1:1000.

Blue: Fluoroshield with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](https://www.genetex.com).