

FPRL1 antibody [HL4245]

Cat. No. GTX643078

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	ICC/IF
Reactivity	Human

Package
100 µl, 25 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
Note : The FcR Blocking Reagent is recommended for use prior to the addition of the researcher's antibody to samples that contain human Fc receptor-expressing cells.	

Not tested in other applications.

Product Note This antibody was raised against human FPRL1 Extracellular domain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide encompassing a sequence within the Extracellular domain of human FPRL1. The exact sequence is proprietary.
Purification	Affinity purified by Protein A.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX643078 ICC/IF Image

FPRL1 antibody [HL4245] detects FPRL1 protein by immunofluorescent analysis.
 Sample: Mock and transfected 293T cells were fixed in ice-cold MeOH for 5 min.
 Green: FPRL1 stained by FPRL1 antibody [HL4245] (GTX643078) diluted at 1:500.
 Blue: Fluoroshield with DAPI (GTX30920).



GTX643078 ICC/IF Image

FPRL1 antibody [HL4245] detects FPRL1 protein by immunofluorescent analysis.
 Sample: HL-60 cells were fixed in 4% paraformaldehyde at RT for 15 min and treated with Fc blocker after blocking step.
 Green: FPRL1 stained by FPRL1 antibody [HL4245] (GTX643078) diluted at 1:100.
 Blue: Fluoroshield™ with DAPI (GTX30920).



For full product information, images and publications, please visit our [website](https://www.genetex.com).